



PRODUCT PROFILE

ET.FOAM

Product	Enzymatic Instrument Pre-Cleaner and Transportation Foam.
Description	ET.Foam is an enzymatic instrument pre-cleaner and transportation foam. ET.Foam prevents blood from drying on instrument surfaces and sustains moisture prior to processing.
Function	ET.Foam effectively cleans away organic matter from surgical instruments. It coats instrument surfaces and penetrates cracks, crevices and lock boxes loosening and dissolving proteinaceous bioburden.
Packaging	1 L spray bottle. Also available is an automated trigger sprayer, which continuously applies thick foam easily.

Features	Benefits
<i>Highly Efficient</i>	On contact ET.Foam begins to loosen and dissolve proteinaceous bioburden.
<i>Neutral pH</i>	ET.Foam with a neutral pH, is a non-corrosive pre-cleaner that effectively dissolves blood, tissue, and other soils.
<i>Time Saving</i>	Instantly begins the enzymatic cleaning process. Perfect for use in surgery, central processing and any department where instrument cleaning is not immediately available. Speeds up cleaning process.
<i>Safe</i>	ET.Foam's bubbly foam prevents spilling and splashing while transporting instruments. Active enzymes begin to break down soils, helping to eliminate hand scrubbing and the risk of personal injury.
<i>User Friendly</i>	Simply spray the ready-to-use solution directly onto soiled instruments. ET.Foam reduces heavy instrument transport weight as compared to liquid-filled trays.

If you have any questions about **ET.Foam**, contact your local area sales representative or call **GERMIPHENE** Corporation toll-free at 1.800.265.9931
Website: www.germiphene.com Email: admin@germiphene.com

