



# PRODUCT PROFILE

## EnVac

<b>Product</b>	Liquid Concentrate Evacuation System Cleaner
<b>Description</b>	An environmentally friendly, highly efficient liquid concentrate for daily cleaning and maintenance of dental aspiration and separation systems.
<b>Function</b>	<b>EnVac</b> will quickly break down organic waste and other residue in the entire evacuation system to keep it thoroughly clean and deodorized.
<b>Packaging</b>	2 Litre Tip 'n' Measure bottle.

Features	Benefits
<i>Easy to Use</i>	Simply dispense 30 ml into 1 litre of water and aspirate through the evacuation system lines daily at the end of the working day.
<i>Effective</i>	<b>EnVac</b> delivers the power to clean your entire system including the traps. It breaks down organic tissue, rope-like saliva and other dental debris where odour-forming bacteria grow.
<i>Cost Effective</i>	A 2 litre container of <b>EnVac</b> provides 66 treatments.
<i>Environmentally Friendly</i>	<b>EnVac</b> is 99.9% biodegradable. Does not contain phenols, aldehydes, phosphates, chlorine or NPE's.
<i>Non-Corrosive</i>	<b>EnVac</b> is compatible with stainless steel, plastic (PVC & ABS), aluminum, and brass. <b>EnVac</b> will not discolour evacuation system materials.
<i>Odour Suppressant</i>	<b>EnVac</b> dissolves odour-causing debris that adheres to the surface of evacuation system lines.

### ***Testimonial from Douglas J.P. Crane, D.D.S.***

*"We had been using an enzyme only vacuum system cleaner, however, when the system was opened for strainer maintenance a very strong sewer gas odour was routinely noticed.*

*Within weeks of using EnVac we noticed large amounts of black biofilm in the strainers and shortly thereafter the system was odour free and free from black biofilm.*

*Our staff found EnVac pleasant, safe and easy to use. I can highly recommend EnVac as a new and well proven product."*

If you have any questions about **EnVac**, contact your local area sales representative or call **GERMIPHENE** Corporation toll-free at 1.800.265.9931  
Website: [www.germiphene.com](http://www.germiphene.com) Email: [admin@germiphene.com](mailto:admin@germiphene.com)

