

# G-NOL

## Disinfectant and Ultrasonic Cleaner

An intermediate-level dual-use solution for cleaning and disinfection of instruments and devices.

G-Nol is a multi-purpose solution which can be used as an instrument disinfectant soak, holding solution, ultrasonic solution, or hard surface disinfectant.

The broad spectrum, intermediate level disinfectant maintains its biocidal activity in the presence of blood and saliva. G-Nol offers complete bactericidal, virucidal, fungicidal and tuberculocidal action in 10 minutes.

G-Nol contains special soaps and surfactants that will thoroughly, safely and quickly clean dental instruments by removing protein and other organic contaminants. The alcohol-free formulation will not discolour or crack office equipment and can be used on metal and rubber instruments.

G-Nol's concentrated formula can be diluted 1:32 with tap water, making it economical to use.

- No activation required
- Complete disinfection in 10 minutes
- Multi-purpose, intermediate-level disinfectant
- Economical, super concentrated
- Pleasantly scented

*A powerful, effective cleaning choice!*



To order this or any Germiphene product, contact us by calling 1-800-265-9931 or via e-mail: [info@germiphene.com](mailto:info@germiphene.com)