



# PRODUCT PROFILE

## COMFORT FIT DISPOSAL GELTRAYS

<b>Product</b>	Disposable Germiphene Geltrays
<b>Description</b>	Basic ingredient: Closed-cell foam. Does not contain natural rubber latex.  Hinged: Dual arch trays for office fluoride treatments Singles: Single arch trays for office fluoride treatments
<b>Function</b>	Proven effective, non-absorbent fluoride tray available in either a dual-arch or single-arch delivery.
<b>Packaging</b>	Hinged: 50 trays available in small, medium or large Singles: 100 trays available in small, medium or large

Features	Benefits
<i>Anatomically Contoured</i>	Anatomically contoured to force fluoride onto all enameled surfaces, particularly interproximal and gingival areas.
<i>Non-Absorbent</i>	Closed-cell foam keeps moisture out and fluoride inside tray. Gel is forced into interproximal and gingival areas.
<i>Effective Design</i>	Distal dam and high walls inhibit gel overflow and fluoride ingestion. Tray depth assures coverage of cervical area. Smooth ridge provides increased patient comfort. Long tray tab ensures clean, easy insertion and removal. Occlusal anatomy forces topical into biting surfaces.
<i>Hinged</i>	Hinged construction enables placement in both arches simultaneously while folding out of the way to reduce gag reflex.
<i>Singles</i>	Upper and lower single trays are dual colour-coded for quick differentiation between sizes.
<i>Latex Free</i>	Ensures no latex allergy reaction.