



# PRODUCT PROFILE

## GLEAM II

<b>Product</b>	Ultrasonic Cleaning Solution
<b>Description</b>	A concentrated, non-ammoniated formulation specifically developed for optimal cleaning efficiency in all types of ultrasonic cleaning units.
<b>Function</b>	<b>GLEAM II</b> is formulated to remove proteinaceous and other organic matter, such as blood, mucus, calculous, etc. from dental instruments.
<b>Packaging</b>	2 litre Tip 'n' Measure bottle and 10 litre container.

Features	Benefits
<i>Easy to use</i>	Simply tip the bottle to fill the upper chamber with the desired amount and pour into ultrasonic tank. This unique dispensing system prevents wastes and spills.
<i>Concentrate</i>	Dilute 1 part of concentrate with 10 parts water.
<i>Superior Cleaning Action</i>	<b>GLEAM II</b> 's formula with powerful detergents and surfactants will clean dental instruments quickly and safely by removing blood and organic soils of all types.
<i>Non-Corrosive</i>	Will not harm stainless steel instruments
<i>Fresh Scent</i>	<b>GLEAM II</b> has a pleasant scent
<i>Reduces Risk Of Cross-Contamination</i>	Ultrasonic cleaning is a critical step in the total infection control process. Prior to sterilization, critical instruments should always be cleaned ultrasonically (to remove all debris and contaminants), thus reducing the risk of cross-contamination associated with manual cleaning.
<i>Protection and Peace of Mind</i>	<b>GLEAM II</b> Ultrasonic Solution plays a vital role in protecting dental professionals by providing superior cleaning results.

If you have any questions about **GLEAM II**, contact your local area sales representative or call **GERMIPHENE** Corporation toll-free at 1.800.265.9931  
Website: [www.germiphene.com](http://www.germiphene.com) E-mail: [admin@germiphene.com](mailto:admin@germiphene.com)

