

GLUTERATE



Sterilant / High Level Disinfectant

Complete TB disinfection in 10 minutes.

Sterilization in 10 hours.



Gluterate has broad-spectrum action, rapidly destroying pathogenic microorganisms. With full contact, Gluterate produces complete inactivation of HIV/AIDS in 30 seconds (with moderate organic soil - 5% blood serum) and viruses in 3 minutes; complete tuberculocidal disinfection in 10 minutes; sterilization in 10 hours.

When used as directed, Gluterate is safe on all medical and dental instruments and devices, lenses, cements, plastics, rubber, steel, aluminum and brass. Gluterate's exclusive anti-corrosive formula prevents rust from developing on metal instruments. Gluterate is the disinfection solution of choice for impressions.

Gluterate requires no activation (unlike alkaline glutaraldehydes) and can be reused for up to 30 days - even in the presence of blood, serum, mucous and soap. Unused portions have an extended shelf life of two years.

- Kills TB in 10 minutes
- Sterilization in 10 hours
- Ready-to-use solution, no activation required
- Clean, fresh mint scent masks the aldehyde odour
- Chemically stable - no degradation at high temperatures
- Reusable up to 30 days, 2-year shelf life on unused portions

Fresh scent, rapid action!



To order this or any Germiphene product, contact us by calling 1-800-265-9931 or via e-mail: info@germiphene.com

When neutralization is necessary - use *Glute Neut* a one-step solution for the disposal of Glutaraldehyde and OPA.

