



PRODUCT PROFILE

HEMASEAL-G **DIN 02246832**

Product	HEMASEAL-G Desensitizing Solution may be applied to all prepared teeth where dentin is involved. Primarily used to desensitize in all restorative situations.
Description	35% Hema (Hydroxyethyl methacrylate), 5% glutaraldehyde. Contains fluoride which is known to reduce the potential for recurrent decay.
Function	Controls dentinal fluid flow by combining with proteins in the tubuli and producing a physiologic barrier.
Packaging	10 mL bottle with dropper dispenser.

Features	Benefits
<i>Versatile</i>	Can be used as a desensitizer alone in cervical abrasions and root sensitivity or in conjunction with all restorative materials including amalgams, crowns, inlays, onlays, veneers, or composite resin restorations. HEMASEAL-G can also be used with any bonding system.
<i>Easy-to-Apply</i>	One-step chairside procedure. Just apply one coat and remove excess with short blast of air. No mixing or curing, no repeat applications required. Provides long-lasting relief in seconds.
<i>Performance</i>	Acts in seconds to seal dentinal tubuli and eliminate fluid shifting, the major cause of post-operative and cervical sensitivity. Use to prevent post-operative and micro-leakage sensitivity. Allows for complete seating of all indirect restorations.
<i>Unique Mode</i>	Glutaraldehyde acts as an anti-microbial, killing most bacteria of Action before sealing the dentinal tubuli. Hema enhances the ability of hydrophilic bonding agents to penetrate into the tubules.
<i>Economical</i>	HEMASEAL-G is only pennies to use per application.
<i>Clinically Proven</i>	Teeth treated with HEMASEAL-G showed a highly significant decrease in sensitivity.

If you have any questions about **HEMASEAL-G**, contact your local area sales representative or call **GERMIPHENE** Corporation toll-free at 1.800.265.9931
Website: www.germiphene.com Email: admin@germiphene.com

