


AN IMPORTANT STEP IN  
PREVENTATIVE MAINTENANCE

got silk?

Instrument  
Protector/  
Lubricant



Scientifically engineered:  
**Silicone Free**  
Corrosion Inhibiting  
Instrument Barrier



# LUBRICATE & PROTECT

## Silk™ Ready-to-Use Spray & Instrument Dip Concentrate

PRE-AUTOCLAVE Instrument Lubricant and Protector aids in the prevention of corrosion of instruments and extends the life of the instruments during the autoclaving process.



NEW	Size	Scent
	4 L                      227 mL	Unscented

**For Use During Autoclaving on:**

- Dental/Surgical Instruments
- Hinges
- Box Locks
- General Equipment

<b>Effective</b>	<i>Bacteriostatic</i>
	<i>Prevents spotting, staining and instrument rust</i>
	<i>Effective in steam or chemical vapour sterilizers</i>
	<i>Routine use reduces cost of instrument repair, replacement &amp; sharpening</i>
	<i>Will <b>not</b> interfere with the microbial killing efficacy of autoclaving</i>
<b>Convenient</b>	<i>Lubricates hinged instruments - restores and maintains movement</i>
	<i>Not sticky, oily or greasy</i>
	<i>Silk™ Ready-to-Use comes in an easy to use spray bottle</i>
<b>Safe</b>	<i>Can be used on a variety of metals: steel, carbide, tungsten, stainless steel and aluminum</i>
	<i>Concentrate requires less storage and less waste in our landfills</i>
	<i>Environmentally Friendly. Biodegradable</i>
	<i>Water based and non-toxic</i>

# GET THE MOST FROM YOUR INSTRUMENTS



For more information contact your District Manager or call Germiphene 1.800.265.9931. Visit us online at [www.Germiphene.com](http://www.Germiphene.com) | Email inquiries: [info@germiphene.com](mailto:info@germiphene.com)

