



PRODUCT PROFILE

SUPER SONIC

Product	Ultrasonic Cleaning Solution
Description	A super concentrated non-ammoniated formulation specifically developed for optimal cleaning efficiency in all types of ultrasonic cleaning units.
Function	SUPER SONIC is formulated to dissolve and remove proteinaceous and other organic matter, such as blood, mucous, and tissue from dental instruments.
Packaging	48 tubes (30 mL each) per box

Features	Benefits
<i>Superior Performance</i>	Provides sparkling cleansing results in less time with less effort. Begins to loosen and dissolve proteinaceous deposits in less than one minute.
<i>Non-Rusting</i>	SUPER SONIC will not harm stainless steel instruments.
<i>Cost-Effective</i>	Requires only 1 tube per 4 L of water. Each box makes 192 L of cleaning solution as needed.
<i>Prevents Build-up</i>	The added chelating agents form a barrier around instruments, protecting them from minerals naturally found in water, which contribute to the build-up of film and deposits.
<i>Bacteriostatic</i>	Prevents bacterial growth.
<i>Versatile</i>	Attractive new packaging is designed with convenience in mind. Boxes are small enough to leave on countertop or easily stored away.
<i>Super Convenient</i>	SUPER SONIC is ready-to-use, no measuring, no waiting, no foaming, no foul odour, and unlike tablets, no residue is left on instruments or in the cleaning tank.

If you have any questions about **SUPER SONIC**, contact your local area sales representative or call **GERMIPHENE** Corporation toll-free at 1.800.265.9931

Website: www.germiphene.com E-mail: admin@germiphene.com

Ultrasonic picture supplied by: Chagrin Safety Supply, Inc., www.chagrinsafetysupply.com.

