

SURGICAL MILK

Concentrated Instrument Lubricant and Rust Inhibitor

Clinically tested to be thoroughly steam penetrable and steam sterilizable.



Surgical Milk is clinically proven to keep instruments in top condition. Routine use reduces the incidence of rust, spots, stains and corrosion on dental and surgical instruments, thus prolonging the life of the instruments.

Surgical Milk pre-autoclave instrument dip will restore movement to hinges as well as ensure proper working of pins and joints. Surgical Milk forms a protective barrier on instruments to help keep them in top working condition.

Surgical Milk is easy to use; no rinsing or wiping is required, simply immerse instruments in the diluted solution - then proceed to sterilize.

- **Bacteriostatic**
- **Effective in steam, dry or chemical vapor sterilizers**
- **Prevents spotting, staining, rusting of instruments**
- **Routine use reduces cost of instrument repair and replacement**
- **Super concentrated – 1 part concentrate to 5 parts water**
- **Does not interfere with the microbial killing efficacy of autoclaving**
- **Water-based – non-oily, non-sticky, non-toxic**

Protect your instruments with Surgical Milk!



To order this or any Germiphene product, contact us by calling 1-800-265-9931 or via e-mail: info@germiphene.com