



PRODUCT PROFILE

GZYME³

Product	Multi Enzymatic Cleaning Solution
Description	GZYME³ is a bacteriostatic triple enzyme cleaner. This three-in-one, non-corrosive cleanser is safe and fast-acting. It can be used as a presoak, ultrasonic detergent and as a linen presoak.
Function	Penetrates and removes organic matter, like blood, proteins, carbohydrates and lipids from endoscopes, medical and dental instruments, minimizing the need to clean instruments by hand.
Packaging	1 L Squeeze'n Measure bottle

GZYME³ is an integral part of the G-Force Infection Control System

Features	Benefits
Highly Efficient	Begins to loosen and dissolve proteinaceous deposits in just one minute.
Superior Performance	GZYME³ uses three enzyme families which are fully effective against carbohydrates, mucopolysaccharides, proteins and lipids.
Environmentally Friendly	Formulated from biodegradable materials, GZYME³ can be disposed via the sink.
Cost Effective	Requires only 10 mL of solution per 4 L of water.
Multiple Application	May be used as an ultrasonic cleaning solution and as an instrument presoak and cleaning solution.
Neutral pH	The GZYME³ low foaming formulation is non-corrosive and compatible with most metals, plastics and rubber.
Powerfully Chelated	The added power of chelating agents form a barrier around instruments, protecting them from minerals naturally found in water that contribute to the build-up of film and binding of hinges.
Bacteriostatic	Prevents bacterial growth.

If you have any questions about **GZYME³**, contact your local area sales representative or call **GERMIPHENE** Corporation toll-free at 1.800.265.9931
 Website: www.germiphene.com E-mail: admin@germiphene.com
 Ultrasonic picture supplied by: Chagrin Safety Supply, Inc., www.chagrinsafetysupply.com.

