

TRÉS CLEAN

Dental Tray and Instrument Cleaner

A mild cleaning formula that leaves trays and instruments sparkling clean.



Unclean

Trés Clean is a concentrated, non-rusting and biodegradable cleaner for use on impression trays, spatulas and other instruments.



Cleaned

Trés Clean will loosen any remaining traces of alginates, zinc-oxide eugenol, carboxylate and phosphate cement and adhesive varnish. The use of warm water accelerates the cleaning action, making residual material easy to rinse or brush off.

The mild cleaning action of Trés Clean is kind to all plastic and non-corroding metal impression trays and instruments.

- Biodegradable
- 1 L makes up to 20 L of use solution
- Safe-to-use on all plastics
- Can be used in ultrasonic machines
- Suited for non-oxidizing metals & plastics

Easily removes the toughest impression material!



To order this or any Germiphene product, contact us by calling 1-800-265-9931 or via e-mail: info@germiphene.com

